

## LIST OF TABLES

<b>Table 3.1.</b> Measurement and Result of 40 % Porous MgO firing temp 1100 °C	28
<b>Table 3.2</b> Measurement and Result of 30% Porous MgO firing temp 1100 °C	28
<b>Table 3.3</b> Measurement and Result of 40 % Porous MgO firing temp 1200 °C	29
<b>Table 3.4</b> Measurement and Result of 30% Porous MgO firing temp 1200 °C	29
<b>Table 3.5</b> Measurement and Result of 40 % Porous MgO firing temp 1300 °C	30
<b>Table 3.6</b> Measurement and Result of 30 % Porous MgO firing temp 1300 °C	30
<b>Table 3.7.</b> Measurement and Result of 40 % Porous MgO firing temp 1100 °C at SP 2h	31
<b>Table 3.8</b> Measurement and Result of 30% Porous MgO firing temp 1100 °C at SP 2h	31
<b>Table 3.9</b> Measurement and Result of 40 % Porous MgO firing temp 1200 °C at soaking period 2 hour.	32
<b>Table 3.10</b> Measurement and Result of 30% Porous MgO firing temp 1200 °C at Soaking period 2 hour	32
<b>Table 3.11</b> Measurement and Result of 40 % Porous MgO firing temp 1300 °C soaking period 2 hour	33
<b>Table 3.12</b> Measurement and Result of 30 % Porous MgO firing temp 1300 °C soaking period 2 hour	33
<b>Table 4.1</b> Measurement and Result of spongy spinel firing temp 1300 °C soaking period 1 h with different mole concentration	42
<b>Table 4.2</b> Measurement and Result of spongy spinel firing temp 1400 °C soaking period 1hour.	46
<b>Table 4.3</b> Measurement and Result of spongy spinel firing temp 1400 °C SP 2h	47
<b>Table 4.4</b> Measurement and Result of spongy spinel firing temp 1300 °C soaking period 2 h	48
<b>Table 5.1</b> Measurement of Relaxation process in different solution treatment with O.1M NaCl	76
<b>Table 5.2</b> Measurement of Relaxation process in different solution treatment with O.2 M NaCl	77
<b>Table 5.3</b> Failure Mode effect analysis template-Minimize Error 6 mol% Solution	82
<b>Table 5.4</b> Failure Mode effect analysis template-Minimize Error for 8 mol% Solution	83
<b>Table 5.5</b> Failure Mode effect analysis template-Minimize Error for 10 mol% Solution	84
<b>Table 9.2.1</b> Corrosion for different bath	120